



STAINLESS STEEL CHILDREN BED BHBD-101

STAINLESS STEEL CHILDREN BED BHBD-101

Children medical Bed built to the most current standards while also meeting the clinician's needs which allows easy care and disinfection. The bed is equipped with Trendelenburg and reverse Trendelenburg positions. Simple and modern design for safety and protection of the child.

Used in Medical unit, Nursing.

Also known as Children Hospital Bed, Medical bed, Hospital bed, Medical child bed, Pediatric bed, Medical pediatric bed, Hospital pediatric bed, Children medical bed, Children pediatric hospital bed.

BHBD-101 STAINLESS STEEL CHILDREN BED

Stainless steel framework
Stainless steel side rail, drop type
Noiseless castors
With mattress
With debris basket

SPECIFICATIONS

Model	BHBD-101
IV pole	Yes
Dimension	750x450x610-1030 mm



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com