



## CHILDREN BED WITH TWO FUNCTION BHBD-107

## CHILDREN BED WITH TWO FUNCTION BHBD-107

Children medical Bed built to the most current standards while also meeting the clinician's needs which allows easy care and disinfection. The bed is equipped with Trendelenburg and reverse Trendelenburg positions. Simple and modern design for safety and protection of the child.

Used in Medical unit, Nursing.

Also known as Children Hospital Bed, Medical bed, Hospital bed, Medical child bed, Pediatric bed, Medical pediatric bed, Hospital pediatric bed, Children medical bed, Children pediatric hospital bed.

## BHBD-107 CHILDREN BED WITH TWO FUNCTION



## SPECIFICATIONS

Model	BHBD-107
Load Bearing	150 kg
Dimension	1770x900 mm
Packing Dimension	wooden carton: 1880x980x1450 mm (assembled)
carton	
1880x1020x310 mm (detached)	Net Weight



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)