



ELECTRIC LIFT PATIENT TRANSFER WHEELCHAIR BHBD-1106

ELECTRIC LIFT PATIENT TRANSFER WHEELCHAIR BHBD-1106

Wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of Controller and other functions reflect the user-friendly design, and provide great convenience for the user.

Used in Patient transfer, patient lift, toilet-bath chair, Patient Indoor-outdoor locomotion.

Also known as Motorized wheelchair, Electric-powered wheelchair.

BHBD-1106 ELECTRIC LIFT PATIENT TRANSFER WHEELCHAIR



Thickened aluminium alloy body, strong bearing capacity

180° opening and closing, easy to realize the transfer of various scenes

Electric motor, one-button button control, easy lifting operation, convenient and quick displacement

Soft and comfortable armrest, comfortable to use

360-degree universal wheel, with locking mechanism, shifting is more secure

PU leather seat cushion, waterproof, can be showered, hollow design, a toilet can be placed below, Easy to extract

SPECIFICATIONS

Model	BHBD-1106
Material	Steel ,PU, Silicone cushion
Maximum loading weight	150 kg
Lifting Range	Seat:30 cm, from 40 cm to70 cm.
Wheel	Two front wheels are with brake
Door width, chair can pass it	At least 58 cm
It suites for bed	Height of bed from 9 cm to 70 cm
Waterproof level	Yes, IP44, top grade waterproof
Control	Remote control
Handles	All foldable
Dimensions	65x57x73 cm
Rated power	96 W
Voltage	DC 24 V
Power supply	Battery, there is a button on battery to power it on /off



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: contact@biolabscientific.com | Website: www.biolabscientific.com