



ELECTRIC TOILET LIFT BHBD-1202

ELECTRIC TOILET LIFT BHBD-1202

Electric Toilet lift is mainly suitable for the elderly, pregnant women, disabled and other people who have inconvenient movements

Used in Nursing home, rehabilitation center, Hospital.

Also known as Toilet lift, Patient transfer chair.

BHBD-1202 ELECTRIC TOILET LIFT



Adjustable Lifting: Adjustable lifting angle for seat height comfort at different heights
Strong Thrust: Double push strength up to 200kg, satisfying 150kg weight and sitting freely

Safe Voltage : DC 24V safe voltage, ultra-low energy consumption

Silent Operation: The whole process of silent operation ensures the comfort of sitting

Self-protection System: When a unilateral failure occurs, the machine automatically starts the self-protection function and stops running

One-button Reset: After use, long press the switch key to reset with one button, while locking other buttons to ensure safe use

SPECIFICATIONS

| | |
|-------------------|--|
| Model | BHBD-1202 |
| Material | Iron, Plastic |
| Seat height | From 45 to 75 cm |
| Weight Capacity | 150 kg |
| Noise | Almost silent electric motor when operation |
| Inclination angle | 15°-16° |
| Size | 57cm width, 65cm length, 47cm height |
| For person | Bariatric person, elderly person, patients disabled person, and pregnant women |
| Rated power | 96 W / 2 A |
| Voltage | DC 24 V |
| Power source | Battery power, rechargeable battery |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com