



ELECTRIC GYNECOLOGICAL TABLE BHBD-210

ELECTRIC GYNECOLOGICAL TABLE BHBD-210

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery. Used in Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation. Also known as Gyneacological table, Gyneacological bed, Examination table, Examination bed, Delivery bed, Gynaecological Delivery Bed.

BHBD-210 ELECTRIC GYNECOLOGICAL TABLE



Adopt high density memory sponge and is once formed, with no seam
Hydraulic lifting by foot pedal, other actions are controlled by pneumatic
Imported Y shaped sealing ring and durable
Stainless steel concealed auxiliary table
The base is made up of high strength ABS
Equipped with foot brake device which is easy for moving the table

SPECIFICATIONS

Model	BHBD-210
Table Length	1300+50 mm
Table width	550+30 mm
Height Adjustment	(700 950)+50 mm
Trendelenburg/rev-trend	$\geq 3^\circ / \geq 10^\circ$
Dimension of auxiliary platform	(350x450)+20 mm
Backboard downward	$\geq 70^\circ$
Standard Accessory	A pair of arm rods, a pair of leg support, a stainless steel sewage basin



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com