



ELECTRIC BED BHBD-401

ELECTRIC BED BHBD-401

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

BHBD-401 ELECTRIC BED



- Central control silent casters: 150 mm
- 5-pieces Aluminium alloy guard rail
- Detachable PP head & foot board
- Surface electrostatic powder coating
- Bed framework made from cold-rolled steel plates

SPECIFICATIONS

Model	BHBD-401
Back-rest lifting angle	0°~75°
Foot-rest lifting angle	0°~45°
Trendelenburg	0°~10°
Reverse trendelenburg	0°~10°
Load Bearing	250 kg
Bed Surface Height	500~1000 mm
Patient surface	900x1960 mm
Drainage bag holders	Yes
IV pole	Yes



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com