



## THREE FUNCTION ELECTRIC BED BHBD-403

## THREE FUNCTION ELECTRIC BED BHBD-403

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

## BHBD-403 THREE FUNCTION ELECTRIC BED



Steel Framework, top class ABS & PP plastic

3 functions controlled by remote control

Castor with central brake system, wheel diameter: 150mm

## SPECIFICATIONS

Model	BHBD-403
Function	3
Back-rest lifting angle	0° - 70°
Knee Gatch	0°-30°
HI - LOW	665mm - 450 mm
Castor	360° swivel castors (Dia. 150 mm), with central brake under the foot board.
Material	Powder coated steel bed frame, top class ABS & PP plastic
Bed board dimension	1940x900 mm
Safe load capacity	< 250 kg
Standard accessory	IV pole
Dimension	2160Lx1120Wx440-655 mm
Net Weight	139 kg
Patient surface	1935Lx900W mm
Height adjustment	440-655 mm
IV pole	Yes



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)