



## AUTOMATIC BED BHBD-420

## AUTOMATIC BED BHBD-420

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

## BHBD-420 AUTOMATIC BED



With the nurse controller and CPR function

Bed board adopts HPL board, x-ray photograph for whole body is available

2 step adjustment at shank section

## SPECIFICATIONS

Model	BHBD-420
Back-rest lifting angle	0°-65°±5°
Foot-rest lifting angle	0°-35°±5°
Bed panel extension	Manual/180mm
Load Bearing	240 kg
Patient surface	1950x910 mm
Dimension	2170x1080x400~750 mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)