



## AUTOMATIC BED BHBD-421

## AUTOMATIC BED BHBD-421

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

## BHBD-421 AUTOMATIC BED



4 pieces European style ABS guardrail, upside down structure

4 luxury central control castors

Bed panel is made of one-step molding cold rolled steel, with excellent performance in anti-slipping and air permeability

1 pair of European style bed head

5 pcs of motors

Driven by imported medical super-silent DC motor, enabling heart displacement position, shock position, CPR position, and examination position

With side controller, overall waterproof level is IPX4

Adopts high precised weighing sensor, LCD display, real-time monitoring patients' weight

Imported German TENTE medical ultra silent castor, central locking system, easy for directional move, universal move and lock

## SPECIFICATIONS

| Model                   | BHBD-421              |
|-------------------------|-----------------------|
| Back-rest lifting angle | ≥65°                  |
| Foot-rest lifting angle | ≥35°                  |
| Trendelenburg           | 12°                   |
| Reverse trendelenburg   | 15°                   |
| Laternal tilt           | 12°                   |
| Height adjustment       | 630-1030 mm           |
| Dimension               | 2240x1110x630/1030 mm |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)