



## TWO CRANK MANUAL BED BHBD-509

## TWO CRANK MANUAL BED BHBD-509

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

## BHBD-509 TWO CRANK MANUAL BED



Steel framework, top class ABS & PP plastic

Collapsible Aluminium alloy protective railing; Detachable PP head & foot board

4 pcs universal castor, dia. 125 mm, Individual brake

A shoe rack under the bed board

## SPECIFICATIONS

Model	BHBD-509
Function	2
Back-rest	0° - 75°
Knee Gatch	0° - 45°
Height (from bed board to ground)	470 mm
Castor	360° swivel castors (Dia. 150 mm), with individual brake.
Material	Powder coated steel bed frame, top class ABS & PP plastic
Bed board dimension	1910x830 mm
Standard accessory	IV pole
Net Weight	65 kg
Patient surface	1910Lx830W mm
Height adjustment	470 mm
IV pole	Yes
Dimension	2080Lx980Wx470 mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)