



## 4 CRANK MANUAL BED BHBD-511

## 4 CRANK MANUAL BED BHBD-511

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

## BHBD-511 4 CRANK MANUAL BED



Steel framework, top class ABS & PP plastic

PP material head / foot board, detachable and easy-clean

4 pcs drop-down PP protective railing, controlled by gas spring

Central brake castors, dia. 150 mm. Brake pedal locate below the foot board

## SPECIFICATIONS

Model	BHBD-511
Function	5
Back-rest	0°-60°
Knee Gatch	0° - 30°
Trendelenburg	0 - 7°
Reverse trendelenburg	0 - 7°
HI - LOW	665 mm - 505 mm
Castor	360° swivel castors (Dia. 150 mm), with central brake under the foot board.
Material	Powder coated steel bed frame, top class ABS & PP plastic
Bed board dimension	1930x990 mm
Standard accessory	IV pole
Net Weight	117 kg
Overall dimension	2120x1055x550-665 mm
Height adjustment	755-485 mm
IV pole	Yes
Dimension	2200Lx1050Wx485-755 mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)