



## SINGLE CRANK HOSPITAL BED BHBD-522

## SINGLE CRANK HOSPITAL BED BHBD-522

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

## BHBD-522 SINGLE CRANK HOSPITAL BED



## SPECIFICATIONS

|                   |                     |
|-------------------|---------------------|
| Model             | BHBD-522            |
| Load Bearing      | 240 kg              |
| Patient surface   | 2150x950 mm         |
| Height adjustment | 500 mm (from floor) |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)