



**MANUAL BED FOUR CRANK 5 FUNCTION  
BHBD-526**

## MANUAL BED FOUR CRANK 5 FUNCTION BHBD-526

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

## BHBD-526 MANUAL BED FOUR CRANK 5 FUNCTION



The main bed and auxiliary bed are connected with each other, auxiliary table is easy for dismantling. The auxiliary table with four wheels

The main bed adjusts various movements for the gears and auxiliary bed is ratched adjustment

Various Accessories are complete

## SPECIFICATIONS

Model	BHBD-526
Table Length	1900+50 mm
Table width	600+20 mm
Height Adjustment	750+50 mm
Back board upwards	≥50°
Backboard downward	≥15°
Seat board downwards	≥8°
Seat board upwards	≥30°
Leg board upwards	≥30°
Leg board downwards	≥30°
Standard accessory	A pair of arm support, a pair of leg support, a pair of arm rods, a pair of foot pedal



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)