



## MANUAL TRANSFER STRETCHER BHBD-815

## MANUAL TRANSFER STRETCHER BHBD-815

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

## BHBD-815 MANUAL TRANSFER STRETCHER



Turn by deasil, the stretcher will be up, and turn by widdershins, the stretcher will be down Bright arrow symbols direct the users how to operate clearly

Central-lock castors, diameter 150mm. Matched with direction castor under the chassis, which help users easy control the stretcher direction

Standard accessories: Mattress, I.V.Pole

## SPECIFICATIONS

Model	BHBD-815
Back-rest	0 - 85°
Height range (bed board to ground)	862 - 566 mm
Standard accessory	Mattress, IV pole
Bed Board	1906x610 mm
Castor	Central lockable 360° swivel castors (Dia. 150 mm)
Overall dimension	2013x700 mm
Net Weight	82 kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)