



CONNECTING TRANSFER STRETCHER BHBD-817

CONNECTING TRANSFER STRETCHER BHBD-817

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

BHBD-817 CONNECTING TRANSFER STRETCHER



It is a manual Connecting Transfer Stretcher, can be used to defend cross-infection

Cart A inside the operation room, Cart B outside, adjusting them to the same height, then push Cart B to Cart A, two carts will be connected, following the orbit, bedboard C can be pushed from one cart to another, and VICE VERSA

Bright arrow symbols direct the users how to operate clearly

Standard accessories: Mattress, I.V. pole

SPECIFICATIONS

Model	BHBD-817
Back-rest	0 - 85°
Height range (bed board to ground)	660 - 910 mm
Standard accessory	Mattress, IV pole
Bed Board	1906x610 mm
Castor	Central lockable 360° swivel castors (Dia. 150 mm)
Overall dimension	3870x840 mm
Net Weight	139 kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com