



## HYDRAULIC TRANSFER STRETCHER BHBD-818

## HYDRAULIC TRANSFER STRETCHER BHBD-818

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

## BHBD-818 HYDRAULIC TRANSFER STRETCHER



It is a type of hydraulic transfer stretcher

Five functions available by pressing the foot-pedals

Right arrow symbols direct the users how to operate clearly

Central-lock castors, diameter 150mm. Matched with direction castor under the chassis, which help users easily control the stretcher direction

Standard accessories: Mattress, I.V.Pole

## SPECIFICATIONS

Model	BHBD-818
Back-rest	0 - 85°
Height range (bed board to ground)	813 - 590 mm
Standard accessory	Mattress, IV pole
Bed Board	1906x610 mm
Castor	Central lockable 360° swivel castors (Dia. 150 mm)
Trendelenburg	0 - 12°
Reverse trendelenburg	0 - 12°
Overall dimension	1906x743 mm
Net Weight	82 kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)