



SLIDE PRINTER BHTP-502

SLIDE PRINTER BHTP-502

Slide printer eliminates the need for writing and labeling by directly printing on the slides with precision. Used in Pathology, Cytology, Histology, Clinical Study Trials, Tracking of specimens. Also known as Slide printer, Cassettes printer.

BHTP-502 SLIDE PRINTER



Laser printing and easy for maintaining. Capability of linking with the management system of the pathology department. Print online or off-line. The writing is clear and delicate, never wears out or fades, which meets the characteristic requirements of keeping slides for long-term reference. No need for preheating, can be immediately used after starting up, safe and noiseless, and free of Maintenance.

SPECIFICATIONS

Model	BHTP-502
Capacity	150 Slides
Average printing speed	$\leq 3S$ / slide
Printing accuracy	≥ 1250 pi
Maximum Capacity	≥ 150 slides
Minimum Capacity	≥ 7 slides
Ambient Humidity	45 - 85 %
Dimensions	465x200x361 mm (LengthxWidthxHeight)
Laser Power	Pulse laser ≤ 2 W
System interface	Integrate with LIS/HIS/P IS
Printing Methods	Print online or off-line
Loading method	Cartridge clip type
External device	Support Bar-code scanner inputting
Consumable items	5×75 mm, 1 - 1.2 mm glass slide



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com