



# BIOLOGICAL INDICATOR INCUBATOR BIBI-102

## BIOLOGICAL INDICATOR INCUBATOR BIBI-102

Biological indicator incubator- Series helps assure the sterility of instruments, vials, media and other products in very short time. It also provides easy evaluation of incubation temperature.

Used in Environmental protection, Research institutes, University laboratories.

Also known as Laboratory Incubator.

## BIBI-102 BIOLOGICAL INDICATOR INCUBATOR

Designed for incubation of EO biological indicators

Customized wells for placement of biological indicator vials

Easy to set incubation temperature

Aluminium blocks available to meet different experimental requirements

## SPECIFICATIONS

Model	BIBI-102
Temperature Range	RT+5-100°C
Temperature Accuracy	0.35°C at 57°C
Temperature Uniformity	0.5°C
Temperature Resolution	0.1°C
Time for Heating	≤6min (20°C to 57°C)
Timer Range	1min-99h 59min
Exterior Dimension	120Wx152Dx112H mm
Display	LCD
Weight	0.85/1.4 kg
Power	120 W
Power Supply	240V, 60Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)