



BULLET ICE MAKER BIBU-101

BULLET ICE MAKER BIBU-101

Biolab Bullet Ice Maker is compact in size, made of highly durable material with energy efficient features. Increased efficiency gives quick production of bullet shaped ice cubes with ice making capacity from 20 kg/24h up to 120 kg/24h. Self contained storage bin eliminates the need of a separate storage unit.

Used in Beverage industries, Biochemical and biological studies, Refrigerated perishables, Food industries, Medical laboratories, Healthcare, Pharmaceutical.

Also known as Ice Maker, Laboratory Bullet Ice Maker.

BIBU-101 BULLET ICE MAKER



- Vertical stirring engine model
- Refrigeration system produces uniform shape ice bullets
- Bell-shaped ice with a hexagon hole
- Highly durable stainless steel 304 2B built-up
- Automatic controller increases working reliability
- Large ice storage bin
- Eco-friendly, energy saving non-Freon refrigerant

SPECIFICATIONS

Model	BIBU-101
Capacity	20 kg/24 h
Storage Bin Capacity	6 kg
Refrigerant	R134a
Condenser Type	Air Cooled
Overall Dimension	330x470x605 mm
Packing Dimension	405x555x665 mm
Weight	28/32 kg
Power	220 W
Power Supply	220V /110 V, 50 / 60 Hz

OPTIONAL ACCESSORIES

Accessory Code	Name
2700606006	Water filter



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com