



## BULLET ICE MAKER BIBU-102

## BULLET ICE MAKER BIBU-102

Biolab Bullet Ice Maker is compact in size, made of highly durable material with energy efficient features. Increased efficiency gives quick production of bullet shaped ice cubes with ice making capacity from 20 kg/24h up to 120 kg/24h. Self contained storage bin eliminates the need of a separate storage unit.

Used in Beverage industries, Biochemical and biological studies, Refrigerated perishables, Food industries, Medical laboratories, Healthcare, Pharmaceutical.

Also known as Ice Maker, Laboratory Bullet Ice Maker.

## BIBU-102 BULLET ICE MAKER



- Vertical stirring engine model
- Refrigeration system produces uniform shape ice bullets
- Bell-shaped ice with a hexagon hole
- Highly durable stainless steel 304 2B built-up
- Automatic controller increases working reliability
- Large ice storage bin
- Eco-friendly, energy saving non-Freon refrigerant

## SPECIFICATIONS

Model	BIBU-102
Capacity	25 kg/24 h
Storage Bin Capacity	7 kg
Refrigerant	R134a
Condenser Type	Air Cooled
Overall Dimension	400x543x670 mm
Packing Dimension	470x615x730 mm
Weight	30/33 kg
Power	240 W
Power Supply	220V /110 V, 50 / 60 Hz

## OPTIONAL ACCESSORIES

Accessory Code	Name
2700606006	Water filter



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)