



## CUBE ICE MAKER BICU-112

## CUBE ICE MAKER BICU-112

Biolab Cube Ice Maker is compact in size, made of highly durable material with energy efficient features. Increased efficiency gives quick production of cube shaped ice cubes with ice making capacity from 20 kg/24h up to 120 kg/24h . Self contained storage bin eliminates the need of a separate storage unit.

Used in Healthcare, Food industries, Refrigerated perishables, Biochemical and biological studies, Medical laboratories, Beverage industries, Pharmaceutical.

Also known as Ice Cuber Machine, Laboratory Cube Ice Maker.

## BICU-112 CUBE ICE MAKER



Micro-computer controlled system with high quality water pump  
Highly durable stainless Steel 201 and 301 construction  
Vertical ice plates  
Highly vigorous ice-making speed with energy efficiency  
Uniform cube shaped ice produced rapidly  
Self contained storage bin

## SPECIFICATIONS

Model	BICU-112
Capacity	300 kg/24 h
Storage Bin Capacity	210 kg
Refrigerant	R134a / R404a
Condenser Type	Water cooled
Machine Head	760x600x560 mm
Cabinet	760x800x980 mm
Packing Machine Head	845x705x850 mm
Packing Cabinet	800x850x1060 mm
Weight	98/107 kg
Power Supply	220 V, 60 Hz

## OPTIONAL ACCESSORIES

Accessory Code	Name
2700606006	Water filter



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)