



## ICE FLAKER BIFL-202

## ICE FLAKER BIFL-202

Biolab Flake Ice maker is compact in size, made of highly durable material with energy efficient features. Increased efficiency gives quick production of ice flakes with ice making capacity from 50 kg/24h up to 130 kg/24h. Self contained storage bin eliminates the need of a separate storage unit.

Used in Pharmaceutical, Refrigerated perishables, Food industries, Biochemical and biological studies, Medical laboratories, Beverage industries, Healthcare.

Also known as Flake Ice Maker, Ice Machine, Laboratory Ice Flaker.

## BIFL-202 ICE FLAKER



- High quality stainless steel construction
- Fully automatic operation system
- Fluorine-free foam thermal cover ensures good insulation
- High speed compact clear flake ice-production
- Advanced self-checking mode function
- Power-overloading, water shortage and ice storage full protection
- Large output and energy efficient
- Easy to clean storage bin

## SPECIFICATIONS

Model	BIFL-202
Capacity	30 kg/24 h
Storage Bin Capacity	10 kg
Refrigerant	R134a
Condenser Type	Air cooled
Overall Dimension	330x470x605 mm
Packing Dimension	405x555x665 mm
Weight	28/31 kg
Power	300 W
Power Supply	220 V, 60 Hz

## OPTIONAL ACCESSORIES

Accessory Code	Name
2700606006	Water filter



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)