



ICE FLAKER BIFL-208

ICE FLAKER BIFL-208

Biolab Flake Ice maker is compact in size, made of highly durable material with energy efficient features. Increased efficiency gives quick production of ice flakes with ice making capacity from 50 kg/24h up to 130 kg/24h. Self contained storage bin eliminates the need of a separate storage unit.

Used in Pharmaceutical, Refrigerated perishables, Food industries, Biochemical and biological studies, Medical laboratories, Beverage industries, Healthcare.

Also known as Flake Ice Maker, Ice Machine, Laboratory Ice Flaker.

BIFL-208 ICE FLAKER



- High quality stainless steel construction
- Fully automatic operation system
- Fluorine-free foam thermal cover ensures good insulation
- High speed compact clear flake ice-production
- Advanced self-checking mode function
- Power-overloading, water shortage and ice storage full protection
- Large output and energy efficient
- Easy to clean storage bin

SPECIFICATIONS

Model	BIFL-208
Capacity	100 kg/24 h
Storage Bin Capacity	40 kg
Refrigerant	R134a
Condenser Type	Air cooled
Overall Dimension	500x611x940 mm
Packing Dimension	570x685x1000 mm
Weight	58/63 kg
Power	520 W
Power Supply	220 V, 60 Hz

OPTIONAL ACCESSORIES

Accessory Code	Name
2700606006	Water filter



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com