



## ICE FLAKER BIFL-209

## ICE FLAKER BIFL-209

Biolab Flake Ice maker is compact in size, made of highly durable material with energy efficient features. Increased efficiency gives quick production of ice flakes with ice making capacity from 50 kg/24h up to 130 kg/24h. Self contained storage bin eliminates the need of a separate storage unit.

Used in Pharmaceutical, Refrigerated perishables, Food industries, Biochemical and biological studies, Medical laboratories, Beverage industries, Healthcare.

Also known as Flake Ice Maker, Ice Machine, Laboratory Ice Flaker.

## BIFL-209 ICE FLAKER



- High quality stainless steel construction
- Fully automatic operation system
- Fluorine-free foam thermal cover ensures good insulation
- High speed compact clear flake ice-production
- Advanced self-checking mode function
- Power-overloading, water shortage and ice storage full protection
- Large output and energy efficient
- Easy to clean storage bin

## SPECIFICATIONS

|                      |                 |
|----------------------|-----------------|
| Model                | BIFL-209        |
| Capacity             | 120 kg/24 h     |
| Storage Bin Capacity | 40 kg           |
| Refrigerant          | R134a           |
| Condenser Type       | Air cooled      |
| Overall Dimension    | 500x611x940 mm  |
| Packing Dimension    | 570x685x1000 mm |
| Weight               | 58/63 kg        |
| Power                | 550 W           |
| Power Supply         | 220 V, 60 Hz    |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)