



GENERAL PURPOSE INCUBATOR BIGP-302

GENERAL PURPOSE INCUBATOR BIGP-302

Biolab Scientific offers a range of high quality, efficient and reliable general purpose incubators that provide superior degree of temperature uniformity. Equipped with microprocessor PID controller, this series of incubators provide exceptional temperature control with added security of over temperature protection. These incubators serve the purpose for most biological analysis, research and general laboratory applications.

Used in Quality Control, Medicine Processing Industries, Research, Pharmaceuticals, Microbial Culturing.
Also known as Lab Incubator, Laboratory General Purpose Incubator.

BIGP-302 GENERAL PURPOSE INCUBATOR



Micro-computer intelligent control system with PID program

Large screen digital display with timer function and touch operation

High quality electrostatic sprayed cold rolling steel exterior and stainless steel inner chamber

Easy to clean semicircular transition chamber with adjustable shelves

Double door sealing for convenient opening

Over temperature protection and auto power shutdown function

SPECIFICATIONS

Model	BIGP-302
Capacity	50 L
Temperature Range	RT +5-65°C
Ambient Temperature	5-40°C
Temperature Accuracy	±0.5°C
Temperature Fluctuation	≤±0.5°C
Temperature Stability	±0.5°C
Temperature Uniformity	±0.5°C
Temperature Resolution	0.1°C
Timer Range	1-9999 min
Load per shelf	15 kg
Internal Dimension	400Wx350Dx350H mm
Exterior Dimension	530Wx480Dx680H mm
Shelves	2
Weight	32/35 kg
Power	250 W
Power Supply	220V, 60Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com