



GENERAL PURPOSE INCUBATOR BIGP-401

GENERAL PURPOSE INCUBATOR BIGP-401

Biolab Scientific offers a range of high quality, efficient and reliable general purpose incubators that provide superior degree of temperature uniformity. Equipped with microprocessor PID controller, this series of incubators provide exceptional temperature control with added security of over temperature protection. These incubators serve the purpose for most biological analysis, research and general laboratory applications.

Used in Quality Control, Medicine Processing Industries, Research, Pharmaceuticals, Microbial Culturing.
Also known as Lab Incubator, Laboratory General Purpose Incubator.

BIGP-401 GENERAL PURPOSE INCUBATOR



- LCD PID control system
- International brand cooling fan
- Stainless steel chamber
- Cavity preheating technology and duct cycling system

SPECIFICATIONS

Model	BIGP-401
Capacity	50 L
Temperature Range	RT +5 ~ 65°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±0.5%
Internal Dimension	380Wx380Dx350H mm
Exterior Dimension	580Wx530Dx600H mm
Shelves	2
Weight	31 kg
Power Supply	220 V 50 Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com