



GENERAL PURPOSE INCUBATOR BIGP-403

GENERAL PURPOSE INCUBATOR BIGP-403

Biolab Scientific offers a range of high quality, efficient and reliable general purpose incubators that provide superior degree of temperature uniformity. Equipped with microprocessor PID controller, this series of incubators provide exceptional temperature control with added security of over temperature protection. These incubators serve the purpose for most biological analysis, research and general laboratory applications.

Used in Quality Control, Medicine Processing Industries, Research, Pharmaceuticals, Microbial Culturing.
Also known as Lab Incubator, Laboratory General Purpose Incubator.

BIGP-403 GENERAL PURPOSE INCUBATOR



- LCD PID control system
- International brand cooling fan
- Stainless steel chamber
- Cavity preheating technology and duct cycling system

SPECIFICATIONS

Model	BIGP-403
Capacity	160 L
Temperature Range	RT +5 ~ 65°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±0.5%
Internal Dimension	600Wx600Dx450H mm
Exterior Dimension	800Wx750Dx700H mm
Shelves	2
Weight	57 kg
Power Supply	220 V 50 Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com