



GENERAL PURPOSE INCUBATOR BIGP-604

GENERAL PURPOSE INCUBATOR BIGP-604

Biolab Scientific offers a range of high quality, efficient and reliable general purpose incubators that provide superior degree of temperature uniformity. Equipped with microprocessor PID controller, this series of incubators provide exceptional temperature control with added security of over temperature protection. These incubators serve the purpose for most biological analysis, research and general laboratory applications.

Used in Quality Control, Medicine Processing Industries, Research, Pharmaceuticals, Microbial Culturing.
Also known as Lab Incubator, Laboratory General Purpose Incubator.

BIGP-604 GENERAL PURPOSE INCUBATOR

Microprocessor controller with Timer function

Natural convection (LGPI-601) and forced convection (Model LGPI-602, LGPI-603, LGPI-604) circulation system maintains temperature stability

High quality stainless steel interior and exterior

Adjustable space between the shelves

Auto-controlled fan speed prevents damage of samples

Self-check function for easy identification of problems

Over temperature and temperature difference alarm

SPECIFICATIONS

Model	BIGP-604
Capacity	90 L
Temperature Range	RT +5-80°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±0.5°C
Temperature Stability	±0.2°C
Temperature Uniformity	±0.5°C
Temperature Resolution	0.1°C
Timer Range	1-9999 min
Internal Dimension	400Wx430Dx520H mm
Exterior Dimension	730Wx680Dx750H mm
Shelves	2
Weight	41/47 kg
Power	600 W
Power Supply	220V, 50Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com