



## ICE LINED REFRIGERATOR BILR-101

## ICE LINED REFRIGERATOR BILR-101

Ice-Lined refrigerators are designed for safe storage of vaccines, pharmaceutical products, blood bag etc. at require temperature 2°C to 8°C. The PFU insulation ensures long hold-over time in case of power failure. These refrigerators are widely applicable in epidemic prevention services, clinics, research institute, hospitals, blood bank, and other allied places.

Used in Pharmaceutical, Chemical, Microbiology, Research Institutions, Immunohematology Department, Blood Bank. Also known as Laboratory Refrigerator.

## BILR-101 ICE LINED REFRIGERATOR

### SPECIFICATIONS

Model	BILR-101
-------	----------



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)