



LABORATORY REFRIGERATOR BLAR-108

LABORATORY REFRIGERATOR BLAR-108

Laboratory refrigerators supplied with digital controller, glass door, alarm system to meet the demanding requirements of laboratory research. These refrigerator has temperature ranging from +2°C to +8°C and all models are combined with two exact sensors and auto defrost.

Used in Industry, Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

Also known as Laboratory Medical Refrigerator.

BLAR-108 LABORATORY REFRIGERATOR

Front opening lockable door with full height handle

Display accuracy to 0.1°C

Audible and Visual alarms: High and low temperature alarm, Sensor failure alarm, Door open alarm, Power failure alarm and low battery alarm.

Interior is made up of high quality ABS material, Exterior is of high quality structural steel, resistant to corrosion and easy to clean

keyboard lock and password protected to prevent the adjustment of the operating parameters

4 units Casters for easy handling

Equipped with back up battery, could support power for alarm and display system for 72 hours during power failure

2-layer glass door with heater to prevent from condensation, especially efficient in the area with high humidity.

Standard with USB interface for temperature data download

SPECIFICATIONS

Model	BLAR-108
Capacity	395 L
Temperature Range	2~8°C
Refrigerant	R600a
Refrigeration System	Huayi Compressor
Defrost Method	Auto Defrost
Controller	Microprocessor Control with LCD Display
Compressor	Huayi Compressor
Insulation	110 mm
Shelves	Adjustable 5 shelves +1 (Storage basket)
Internal Size	550x540x1415 mm
External Size	650x697x1815 mm
Weight	140 / 146 Kg
Power Supply	220V,50/ 60HZ, 110V,50/60HZ



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com