



# LABORATORY REFRIGERATOR BLAR-202

## LABORATORY REFRIGERATOR BLAR-202

Laboratory refrigerators supplied with digital controller, glass door, alarm system to meet the demanding requirements of laboratory research. These refrigerator has temperature ranging from +2°C to +8°C and all models are combined with two exact sensors and auto defrost.

Used in Industry, Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

Also known as Laboratory Medical Refrigerator.

## BLAR-202 LABORATORY REFRIGERATOR

Safety door lock

Auto defrost function

High ambient temperature and humidity area

Wide voltage from 192 V to 242 V

No condensation glass door

Functions: High / Low Temperature, High / Low Temperature Recorder, Power Failure, Low Battery, Door Ajar, Remote Alarm

Optional Accessories: Foaming door, Temperature Recorder

## SPECIFICATIONS

Model	BLAR-202
Capacity	110 L
Temperature Range	2~8°C
Cabinet Type	Under-counter
Cooling Performance	4°C
Temperature Resolution	0.1°C
Refrigerant	R134a, CFC free ( R600a optional )
Defrost Method	Auto Defrost
Cooling Type	Forced air cooling
Material	Galvanized steel powder coating
Shelves	2
Climatic class	ST
Internal Size	516Wx388Dx585H mm
External Size	600Wx531Dx805H mm
Power	80 W
Power Supply	220V,50/ 60HZ, 110V,50/60HZ



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)