



## LABORATORY REFRIGERATOR BLAR-205

## LABORATORY REFRIGERATOR BLAR-205

Laboratory refrigerators supplied with digital controller, glass door, alarm system to meet the demanding requirements of laboratory research. These refrigerator has temperature ranging from +2°C to +8°C and all models are combined with two exact sensors and auto defrost.

Used in Industry, Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

Also known as Laboratory Medical Refrigerator.

## BLAR-205 LABORATORY REFRIGERATOR

Safety door lock

Auto defrost function

High ambient temperature and humidity area

Wide voltage from 192 V to 242 V

No condensation glass door

Functions: High / Low Temperature, High / Low Temperature Recorder, Power Failure, Low Battery, Door Ajar, Remote Alarm

Optional Accessories: Foaming door, Temperature Recorder

## SPECIFICATIONS

| Model                  | BLAR-205                           |
|------------------------|------------------------------------|
| Capacity               | 380 L                              |
| Temperature Range      | 2~8°C                              |
| Cabinet Type           | Vertical                           |
| Cooling Performance    | 4°C                                |
| Temperature Resolution | 0.1°C                              |
| Refrigerant            | R134a, CFC free ( R600a optional ) |
| Defrost Method         | Auto Defrost                       |
| Cooling Type           | Forced air cooling                 |
| Material               | Galvanized steel powder coating    |
| Shelves                | 5                                  |
| Climatic class         | ST                                 |
| Internal Size          | 516Wx470Dx1640H mm                 |
| External Size          | 600Wx583Dx1860H mm                 |
| Power                  | 135 W                              |
| Power Supply           | 220V,50/ 60HZ, 110V,50/60HZ        |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)