



LABORATORY REFRIGERATOR BLAR-205

LABORATORY REFRIGERATOR BLAR-205

Laboratory refrigerators supplied with digital controller, glass door, alarm system to meet the demanding requirements of laboratory research. These refrigerator has temperature ranging from +2°C to +8°C and all models are combined with two exact sensors and auto defrost.

Used in Industry, Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

Also known as Laboratory Medical Refrigerator.

BLAR-205 LABORATORY REFRIGERATOR

Safety door lock

Auto defrost function

High ambient temperature and humidity area

Wide voltage from 192 V to 242 V

No condensation glass door

Functions: High / Low Temperature, High / Low Temperature Recorder, Power Failure, Low Battery, Door Ajar, Remote Alarm

Optional Accessories: Foaming door, Temperature Recorder

SPECIFICATIONS

Model	BLAR-205
Capacity	380 L
Temperature Range	2~8°C
Cabinet Type	Vertical
Cooling Performance	4°C
Temperature Resolution	0.1°C
Refrigerant	R134a, CFC free (R600a optional)
Defrost Method	Auto Defrost
Cooling Type	Forced air cooling
Material	Galvanized steel powder coating
Shelves	5
Climatic class	ST
Internal Size	516Wx470Dx1640H mm
External Size	600Wx583Dx1860H mm
Power	135 W
Power Supply	220V,50/ 60HZ, 110V,50/60HZ



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com