



LABORATORY REFRIGERATOR BLAR-411

LABORATORY REFRIGERATOR BLAR-411

Laboratory refrigerators supplied with digital controller, glass door, alarm system to meet the demanding requirements of laboratory research. These refrigerator has temperature ranging from +2°C to +8°C and all models are combined with two exact sensors and auto defrost.

Used in Industry, Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

Also known as Laboratory Medical Refrigerator.

BLAR-411 LABORATORY REFRIGERATOR



SPECIFICATIONS

| | |
|------------------------------|---|
| Model | BLAR-411 |
| Capacity | 316 L |
| Temperature Range | 2~8°C |
| Ambient Temperature | 10~32°C |
| Refrigerant | R290 |
| Defrost Method | Auto Defrost |
| Controller | Microprocessor Control with digital Display |
| Control Sensor | NTC |
| Cooling Type | Forced air cooling |
| Shelf/drawer | Shelves/5 |
| Alarms | High/low temperature, Power failure, Thermostat failure, Door ajar, Power failure backup system(alarm) 8h |
| Insulation | PURF |
| Material Inside | Color sprayed steel |
| Material Outside | Color sprayed steel |
| Noise Level | 50 dB |
| Internal Size | 540x425x1380 mm |
| External Size | 660x555x1915 mm |
| Package Size | 720x680x2040 mm |
| CBM | 0.9 |
| 20GP/40GP/40HQ | 27/54/54 |
| NT./GT. (kg) | 92.5/100kg |
| Power consumption (kW.h/24h) | 1.78 |
| Power | 228W |
| Current (A) | 1.2A |

| | |
|--------------------------|--|
| Voltage/frequency (V/Hz) | 220/50Hz |
| Standard Accessories | Foot, Castor, Glass Low-E type Outer door, Door Lock, Test hole 25mm diameter, LED Lamp, Remote alarm port |
| Optional Accessory | Temperature recorder Datalogger, USB port, RS232/485 port |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com