



DESKTOP LAMINAR AIRFLOW BLDK-102

DESKTOP LAMINAR AIRFLOW BLDK-102

Compact desktop design required for operations in an ultra clean, dust free environment. Small size saves precious laboratory space. Vertical and horizontal laminar flow cabinets available suit your needs.

Used in Medical, Pharmaceutical, Research, Manufacturing areas, Electronic setors, Industrial sectors, Laboratory, Healthcare.

Also known as Benchtop Laminar Air Flow, Laminar Flow Cabinet, Clean Bench, Laboratory Laminar Flow, Laboratory Laminar Air Flow, Tissue Culture Hood, Laminar Air Flow Cabinet, Minimal-Turbulence Air Flow, Laminar Flow Cabinet, Laminar Flow Hood, Laminar Flow.

BLDK-102 DESKTOP LAMINAR AIRFLOW

Air Cleanliness: Class 100 @ $\geq 0.5 \mu\text{m}$
Average Air Flow Velocity: 0.3-0.45 m/s
Fluorescent Lamp: 8 W
Noise Level: $\leq 60 \text{ dB(A)}$

SPECIFICATIONS

Model	BLDK-102
Type	Vertical Laminar Flow
Air Cleanliness	Class 100 @ $\geq 0.5 \mu\text{m}$
Average Air Flow Velocity	0.3-0.45 m/s
Fluorescent Lamp	8 W
Internal Size	600Wx400Dx535H mm
Overall Dimension	650x700x625 mm
Noise Level	$\leq 60 \text{ dB(A)}$
Power	400 W
Power Supply	220 V, 50 Hz

ACCESSORIES

Accessory Code	Name	Description	Unit
2800606007	Power cord	AC220V 60Hz	1
2800606008	Fuse - 5A		5
2800607006	Veritical Clean Bench		1
2800606009	Fluorescent lamp- 14 W		1
2800606010	Ultraviolet lamp - 14 W		1



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com