



## HOMOGENIZER BLHM-2001

# HOMOGENIZER BLHM-2001

## BLHM-2001 HOMOGENIZER



- Designed smart and light
- Can operate in one hand
- Speed adjustable, digital display speed. Speed precisely controlled to avoid splashing
- Suitable for homogeneous dispersion of a small amount of sample in small container
- High quality stainless steel dispersing head, anti-corrosion, high temperature resistance for disinfecting
- Easy to take off the dispersing head for cleaning
- High speed and high efficiency motor, long life
- Optional H-type stand, height and angle can be adjusted, convenient to take and store

## SPECIFICATIONS

Model	BLHM-2001
Speed Range	8 krpm ~ 35 Krpm(1 krpm Increment)
Rotor Line Speed	6.3 m/ sec~14 m/ sec
Working Noise	72 dB
Voltage	AC 220 V, 50 Hz
Dimension	W.46 x D.55 x H.230 mm
Net Weight	0.8 kg
Power Supply	160 W/110 W



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)