



AUTOMATIC SEALING AND CAPPING MACHINE BLIH-103

AUTOMATIC SEALING AND CAPPING MACHINE BLIH-103

Automatic Sealing and Capping Machine is specifically designed to speed up the process of packaging of your product but also ensure the sterility of the product and that the product does not get contaminated
Used in Pharmaceuticals, Mobile inspection laboratories, Disease prevention and control centers, Food Industry, Beverage Industry.
Also known as Screw on capping machine.

BLIH-103 AUTOMATIC SEALING AND CAPPING MACHINE



- Fully automatic and one key operation, no need of additional manual operation
- Adapts to varieties of plate types
- Adjustable sealing pressure
- A variety of communication interfaces are reserved, which can be quickly connected to other machines
- Each roll of film can seal more than 350 plates(reference quantity)

SPECIFICATIONS

Model	BLIH-103
Positions	1pc
Method	Cold Seal
Speed	<16s/plate
Operating Mode	Manual plate change, automatic sealing
Application	PCR plate, deep well plate and micro-plate conforming to SBS standard(Customizable)
Control Mode	Buttons
Interface	USB
External Size	310x380x410 mm
Package Size	454x382x464 mm
Weight	Net Weight- 16.5kg Gross Weight- 20.5kg
75 W	Power Supply



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com