



## SAMPLE TUBE SCANNER BLIH-201

## SAMPLE TUBE SCANNER BLIH-201

Sample tube scanner is specifically designed to operate on top of the sample tubes without applying any weight, avoiding any damage.

Used in Laboratory, Research, Medical.

Also known as Tube Reader, Tube Scanner, Tube Barcode scanner, Sample Scanner.

## BLIH-201 SAMPLE TUBE SCANNER



It can be finish within 5 seconds for one of 96 holes rack or 384 holes rack.

Advanced decoder technology can be identified stable and secure even if the code is fuzzy.

Optional item for 1-D code rack scanning function.

It can be scanned difference of height sample tube in the same rack .

Compact, ergonomic design assures efficient of workflow.

Detection of rack type automatically for setting optimum scanning.

## SPECIFICATIONS

|                     |   |
|---------------------|---|
| Model               | BLIH-201  |
| Rack type available | 13,24,48,96,384 holes rack  |
| 1-D Code available  | 2/5 Industrial / Interleaved, Code 39,Code 128, Pharmacode, Codabar, EAN 13 |
| 2-D Code available  | Datamatrix ECC 200, PDF417, QR Code   |
| Camera              | 10 mega-pixel camera CMOS   |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)