



## PORTABLE ULTRASONIC PROCESSOR BLIH-401

## PORTABLE ULTRASONIC PROCESSOR BLIH-401

Portable Ultrasonic Processor is used to emulsify, disperse, homogenization, extraction, anti-foaming, cleaning, preparation of nano-materials, dispersion, and speed up chemical reactions

Used in Biotechnology, Pharmaceutical, Chemical, Industrial, Environmental.

Also known as Portable Ultrasonic Homogenizer Sonicator, Portable ultrasonic cell crusher.

## BLIH-401 PORTABLE ULTRASONIC PROCESSOR



## SPECIFICATIONS

Model	BLIH-401
Operating frequency range	20-25 KHz, frequency automatic tracking
The end of the horn diameter	Φ 3 or 5
Ultrasonic Power	150W, can be adjusted, hit into ways to fine-tune each of 5W
Parameters sample handling capacity	0.1-100ml (need to select the appropriate probe)
Generator output power	150W (max)
Operating voltage	220V±5%



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)