



Product Image Coming Soon

MULTI-PARAMETER ANALYZER BMET-1104

MULTI-PARAMETER ANALYZER BMET-1104

Multi-parameter Analyzer simultaneously measures Ion, Conductivity, Dissolved Oxygen, Temperature that otherwise requires four different meters. With the appropriate probes installed, it can also simultaneously measure several or one parameters.

BMET-1104 MULTI-PARAMETER ANALYZER



Product Image Coming Soon

SPECIFICATIONS

Model	BMET-1104
DC Voltage	
2 V/20 V/200 V	$\pm (0.5\%+3)$
600 V	$\pm (0.8\%+5)$
AC Voltage	
2 V/20 V/200 V	$\pm (0.8\%+3)$
600 V	$\pm (1.0\%+5)$
Resistance	
2 k Ω /20 k Ω /200 k Ω /2 M Ω /10 M Ω	+ (1.0%+3)
Display	2000 Counts
Data Hold	Yes
Continuity	Yes
NCV	Yes
Work Light	Yes
Power Supply	2x 1.5 V AAA Batteries
Weight	Approx. 133.5/146.5 g
Size	130x62x27 mm
Safety Standard/Rating	EN61010-1,EN61010-2-033,EN61326,CAT.III 600 V



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com