



Product Image Coming Soon

MULTI-PARAMETER ANALYZER BMET-1105

MULTI-PARAMETER ANALYZER BMET-1105

Multi-parameter Analyzer simultaneously measures Ion, Conductivity, Dissolved Oxygen, Temperature that otherwise requires four different meters. With the appropriate probes installed, it can also simultaneously measure several or one parameters.

BMET-1105 MULTI-PARAMETER ANALYZER



Product Image Coming Soon

SPECIFICATIONS

Model	BMET-1105
DC Voltage	
200 mV/2 V/20 V/200 V	± (0.5%+3)
600 V	± (0.8%+5)
AC Voltage	
200 V/600 V	± (0.8%+3)
600 V	± (1.0%+5)
DC Current	
200 µA/2000 µA/20 mA/200 mA	± (1.8%+5)
AC Current	
200 uA/2000 uA/20 mA/200 mA	± (2.0%+5)
Resistance	
200 Ω/2 kΩ/20 kΩ/200 kΩ/2 MΩ/10 MΩ	± (1.0%+3)
20 MΩ	± (1.2%+15)
Display	2000 Counts
Data Hold	Yes
Continuity	Yes
NCV	Yes
Work Light	Yes
Power Supply	2x 1.5 V AAA Batteries
Weight	Approx. 133.5/146.5 g
Size	130x62x27 mm
Safety Standard/Rating	EN61010-1,EN61010-2-033,EN61326,CAT.III 600 V



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com