



MAGNETIC STIRRER BMGS-103

MAGNETIC STIRRER BMGS-103

Magnetic stirrer series offers many work modes, Strong magnetism, 30° bevel control panel, Brushless DC motor for convenient operations. Automatic reverse rotation for better-mixing result. Slip-proof and stable. Programmable memory can store nine common set S1 ~ S9, Meet the different kinds of experimental demands.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Multistation Magnetic Stirrer, Digital Magnetic Stirrer.

BMGS-103 MAGNETIC STIRRER

4-channels

Specially designed for cell culture that needs low speed stirring

Programmable, make the operation mode automatically

Digital LED displays the actual speed

Many work modes, including the constant speed stirring, cycle stirring, two speed stirring and combination cycle and two speed stirring

SPECIFICATIONS

| | |
|--------------------|-------------------|
| Model | BMGS-103 |
| Capacity | 3 L |
| Stirring Positions | 4 |
| Speed Range | 5-200 rpm |
| Dimension | 395x450x90 mm |
| Power | 15 W |
| Power Supply | 100-240V, 50/60Hz |

OPTIONAL ACCESSORIES

| Accessory Code | Name | Description | Capacity |
|----------------|---------------|---------------------------------|----------|
| 3400807006 | Spinner Flask | Standard; Flat Bottom; Complete | 500 ml |
| 3400807007 | Spinner Flask | Standard; Flat Bottom; Complete | 3000 ml |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com