



MAGNETIC STIRRER BMGS-103

MAGNETIC STIRRER BMGS-103

Magnetic stirrer series offers many work modes, Strong magnetism, 30° bevel control panel, Brushless DC motor for convenient operations. Automatic reverse rotation for better-mixing result. Slip-proof and stable. Programmable memory can store nine common set S1 ~ S9, Meet the different kinds of experimental demands.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Multistation Magnetic Stirrer, Digital Magnetic Stirrer.

BMGS-103 MAGNETIC STIRRER

4-channels

Specially designed for cell culture that needs low speed stirring

Programmable, make the operation mode automatically

Digital LED displays the actual speed

Many work modes, including the constant speed stirring, cycle stirring, two speed stirring and combination cycle and two speed stirring

SPECIFICATIONS

Model	BMGS-103
Capacity	3 L
Stirring Positions	4
Speed Range	5-200 rpm
Dimension	395x450x90 mm
Power	15 W
Power Supply	100-240V, 50/60Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	Capacity
3400807006	Spinner Flask	Standard; Flat Bottom; Complete	500 ml
3400807007	Spinner Flask	Standard; Flat Bottom; Complete	3000 ml



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com