



## MAGNETIC STIRRER BMGS-108

## MAGNETIC STIRRER BMGS-108

Magnetic stirrer series offers many work modes, Strong magnetism, 30° bevel control panel, Brushless DC motor for convenient operations. Automatic reverse rotation for better-mixing result. Slip-proof and stable. Programmable memory can store nine common set S1 ~ S9, Meet the different kinds of experimental demands.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Stirrer, Digital Magnetic Stirrer.

## BMGS-108 MAGNETIC STIRRER



Speed range of 0-1500 rpm

Maximum stirring quantity 10 L

Glass ceramic plate, excellent chemical resistance and easy to clean

## SPECIFICATIONS

Model	BMGS-108
Capacity	10 L
Speed Range	0-1500 rpm
Temperature	5-40°C



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)