



## STEREO ZOOM MICROSCOPE BMIC-903

## STEREO ZOOM MICROSCOPE BMIC-903

Stereo Zoom Microscope also known as the dissecting microscope is type of optical microscope that allows the user to see a three-dimensional view of a specimen.

Used in Pharmaceutical, University, Medical, Laboratory, Research.

## BMIC-903 STEREO ZOOM MICROSCOPE



Excellent optical function with Natural light binocular illumination. Innovative stand structure and sharp image display. Convenient operation and special design for viewing incubating cell tissue. Valid mould proof design. It can be used under the high temperature and high humidity environment.

## SPECIFICATIONS

Model	BMIC-903
Illumination	Natural light Binocular
View Head And Body	Straight Head
Objective	2X
Working Distance	80 mm
Eyepiece	WF10X/20
Power Supply	NA
Package Size	300x200x150 mm
Gross Weight	3 kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)