



STEREO ZOOM MICROSCOPE BMIC-904

STEREO ZOOM MICROSCOPE BMIC-904

Stereo Zoom Microscope also known as the dissecting microscope is type of optical microscope that allows the user to see a three-dimensional view of a specimen.

Used in Pharmaceutical, University, Medical, Laboratory, Research.

BMIC-904 STEREO ZOOM MICROSCOPE



Excellent optical function with Tranilluminator / inclined light illumination. Innovative stand structure and sharp image display. Convenient operation and special design for viewing incubating cell tissue. Valid mould proof design. It can be used under the high temperature and high humidity environment.

SPECIFICATIONS

Model	BMIC-904
Illumination	Tranilluminator / inclined light
View Head And Body	Inclined viewing Up-rotatable body
Objective	2X&4X, 90°rotated
Working Distance	57 mm
Eyepiece	WF10X/20
Power Supply	AC110/220 V \pm 10%, 50/60 Hz
Package Size	360x170x290 mm
Gross Weight	6 kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com