



## MULTISTATION STIRRER BMSS-102

## MULTISTATION STIRRER BMSS-102

Powerful Multi-Station magnetic stirrers with 4 to 16 stirring positions to minimize bench space. Uniform agitation reaction conditions are achieved by synchronized operation. Speed ranges from 0-1100 rpm. Individually controlled stirring for consistency with various samples.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Multistation Magnetic Stirrer, Digital Magnetic Stirrer.

## BMSS-102 MULTISTATION STIRRER



Digital display stirring speed at real-time, stir bar operate synchronously and provide uniform agitation reaction conditions.

Efficient Multistation magnetic stirrer configure 4/8/12 stirring positions to use lab space effective.

Low speed but stable, high speed but powerful. Low noise and maintenance free.

Stainless steel is coated by silicone resin film disk.

Non-slip and anti-corrosion.

DC brushless motor.

## SPECIFICATIONS

Model	BMSS-102
Capacity	400 ml
Stirring Positions	8
Speed Range	0~1100 rpm
Platform Material	Stainless steel with silicon resin film
Speed Display	Digital
Motor Type	Brushless DC motor
Dimension (WxDxH)	205x480x50 mm
Fuse	250V,1A, $\phi$ 5x20
Weight	4 kg
Power	20 W
Power Supply	AC100-230V, 50/60Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)