



## NATURAL CONVECTION OVEN BONC-102

## NATURAL CONVECTION OVEN BONC-102

With easy to clean interiors and natural convection circulation system for efficient heating of chamber, this oven can be relied on for general laboratory applications.

Used in Drying, Roasting, Laboratory, Research, Engineering, Industry, Agricultural, Biology.

Also known as Laboratory Oven, Laboratory natural Convection Oven.

## BONC-102 NATURAL CONVECTION OVEN



304 stainless steel, mirror polishing processing, easy to clean and maintain

Natural convection with low noise

PID controller with over temperature alarm and Timer function ensures precise and reliable control

## SPECIFICATIONS

Model	BONC-102
Capacity	56 L
Temperature Range	RT+10-200°C
Ambient Temperature	+5-35°C
Display Resolution	0.1°C
Internal Dimension	400Wx330Dx415H mm
Exterior Dimension	540Wx550Dx720H mm
Shelves	2
Power consumption	1000 W
Electrical Requirement	AC 220V, 50Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)