



## NATURAL CONVECTION OVEN BONC-103

## NATURAL CONVECTION OVEN BONC-103

With easy to clean interiors and natural convection circulation system for efficient heating of chamber, this oven can be relied on for general laboratory applications.

Used in Drying, Roasting, Laboratory, Research, Engineering, Industry, Agricultural, Biology.

Also known as Laboratory Oven, Laboratory natural Convection Oven.

## BONC-103 NATURAL CONVECTION OVEN



304 stainless steel, mirror polishing processing, easy to clean and maintain

Natural convection with low noise

PID controller with over temperature alarm and Timer function ensures precise and reliable control

## SPECIFICATIONS

Model	BONC-103
Capacity	96 L
Temperature Range	RT+10-200°C
Ambient Temperature	+5-35°C
Display Resolution	0.1°C
Internal Dimension	450Wx430Dx505H mm
Exterior Dimension	590Wx650Dx810H mm
Shelves	2
Power Consumption	1400 W
Electrical Requirement	AC 220V, 50Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)