



NATURAL CONVECTION OVEN BONC-104

NATURAL CONVECTION OVEN BONC-104

With easy to clean interiors and natural convection circulation system for efficient heating of chamber, this oven can be relied on for general laboratory applications.

Used in Drying, Roasting, Laboratory, Research, Engineering, Industry, Agricultural, Biology.

Also known as Laboratory Oven, Laboratory natural Convection Oven.

BONC-104 NATURAL CONVECTION OVEN



304 stainless steel, mirror polishing processing, easy to clean and maintain

Natural convection with low noise

PID controller with over temperature alarm and Timer function ensures precise and reliable control

SPECIFICATIONS

| | |
|------------------------|-------------------|
| Model | BONC-104 |
| Capacity | 140 L |
| Temperature Range | RT+10-200°C |
| Ambient Temperature | +5-35°C |
| Temperature Stability | ±1°C |
| Display Resolution | 0.1°C |
| Internal Dimension | 520Wx500Dx575H mm |
| Exterior Dimension | 660Wx720Dx880H mm |
| Shelves | 3 |
| Power Consumption | 2000 W |
| Electrical Requirement | AC 220V, 50Hz |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com